FLOCARE® PUR TUBE with guidewire

Manufacturer: Nutricia Medical Devices BV, Taurusavenue 167, 2132 LS Hoofddorp, The Netherlands, www.nutriciaflocare.com
Production location: Nutricia Pharmaceutical Wuxi Co., Ltd., P.R. China

CE
0344

STERILE EO MD

200593_A
Issue Date
October 2022

FLOCARE® PUR TUBE with guidewire

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
outer cover right page to be folded inside (3)

Supplies needed / Erforderliches zubehör zum legen einer transnasalen
 magensonde / Benodigdheden / Matériel nécessaire à la pose / Udstyr / Utrustning /
 Tarvikkeet / Utstyr / Suministros necesarios / Materiali necessari / Material
 necessário / Potřebné vybavení / Απαιτούμενος εξοπλισμός / Gerekenler /
 Nepieciešamie priekšmeti

The F
 (maxi

The tu
 and d

	Re - F chi - C - A
	Ma
	Col con roto
	X-r
	Gu
	De fluid
	Lat
	Tip

EN • FLOCARE® PUR TUBE with guidewire

PRODUCT DESCRIPTION

The Flocare® Pur nasogastric feeding tube is a transnasal tube recommended for patients who require short term (maximum 6 weeks) administration of enteral feeding directly into the stomach.

The tubes are available in a wide range of lengths and charrières, made of comfortable and skin-friendly materials and designed for an optimal feed flow. The tubes have the following characteristics:

Materials & key features	Flocare® PUR range
Recommended size: - Prematurely & newborn children (neonates) - Children - Adults	: maximum CH6 : CH5-8; : minimum CH 10 <i>The selection of the CH-size is to be decided by the trained healthcare professional while considering patient age, condition, type of feeding and feeding regime.</i>
Materials	Transparent polyurethane tubing, latex and DEHP-free materials
Color coded ENFIT™ connectors equipped with a rotating swivel:	CH5 grey, CH6 light green, CH8 blue, CH10 Black, CH12 white, CH14 green
X-ray visibility	3 white opaque lines
Guidewire	Lubricated guidewire ends before the first lateral opening of the tube leaving the tube tip bend slightly, following the curvature of the pharynx which eases insertion
Dedicated port for gastric fluid drainage	CH 14 variant available
Lateral openings	2 Suction tube: 4
Tip	Open: CH6, CH8, CH10, CH12, CH14 Closed: CH5

INDICATIONS FOR ENTERAL TUBE NUTRITION

Enteral tube nutrition or tube feeding is indicated for patients at nutritional risk or malnourished who:

- Cannot meet their nutrient requirements by normal dietary intake and;
- Have a functioning gastrointestinal tract

INDICATIONS FOR NASOGASTRIC FEEDING TUBES

Indication for therapeutic reasons:

- Intra-gastric enteral feeding
 - Stomach decompression or suction.
- Dedicated tubes with suction ports should be used.

CONTRAINDICATIONS FOR NASOGASTRIC FEEDING TUBES

Nasogastric feeding is unsuitable for patients with severely compromised gastric emptying owing to partial obstruction or ileus. Relative contraindications for nasogastric tubes are comparable to general contraindications for enteral feeding, such as acute pancreatitis, acute abdomen, motility disorders, severe malabsorption and severe Crohn's disease. In these cases, it depends on the judgement of the physician and the clinical situation of the patient whether enteral feeding is given.

Notice: any serious incident that has occurred in relation to this device is to be reported to Nutricia and the local competent authority.


INSTRUCTIONS FOR USE

Supplies needed (see page 4):

- | | |
|--|--|
| 1. Flocare® nasogastric tube | 8. Hand sanitizer |
| 2. Hypoallergenic tape and optional marker | 9. Alcohol |
| 3. 60ml enteral feeding syringe | 10. Paper handkerchief |
| 4. Personal protective equipment (apron, gown) | 11. Protective towel to put equipment in the table |
| 5. Cup of water with straw or ice chips | 12. Protective towel for patient |
| 6. Emesis bowl | 13. Small container of water |
| 7. pH indicator strips | 14. Disposable latex-free gloves |


PREPARATION

△ **Nasogastric tubes should only be inserted by trained healthcare professionals because of the risk of misplacement and oesophageal or pulmonary perforation.**


1. Preparation.

Explain procedure to the patient and gather the supplies needed. Place the patient in a comfortable sitting or semi-recumbent position. If unable, lay the patient on one side.


2. Preparation.

Agree on a signal with the patient that he or she can use to stop the procedure (such as raising the hand).


3. Preparation

Wash and disinfect your hands and put on gloves. Please follow local policies for PPE (Personal Protective Equipment).

4. Tube measurement


Determine the tube length to be inserted; Use the end of the NG-tube to measure the distance from the tip of the nose (**N**), to the ear (**E**), and then to the Xiphoid process of the sternum (**X**). Mark the tube at this point with the marker or tape. Local guidelines can differ in terms of measurement method. Always consult your local guidelines.

5. Guidewire insertion


Fully insert the guidewire in the tube and rotate the swivel to connect to the connector. Do not lubricate the guidewire by flushing water down the tube.

TUBE PLACEMENT


6. Tip lubrication

Lubricate the tip with water, to facilitate insertion


7. Placing the tube

Inspect the patient's nostrils for a deviated septum or visible nasal polyps; use the larger nostril for your insertion. Pass the tube along the bottom of the nostril. Direct the tube posteriorly, aiming tip parallel to nasal septum following the surface of the hard palate. Advance the tube to nasopharynx, allowing tip to seek its own passage.

8. Placing the tube

As the patient swallows sips of cold water (Infants can suck on a pacifier during the procedure), advance the tube through esophagus into stomach with gentle motion.

9. Placing the tube

Push the tube gently, never force it. Continue until the measured mark is reached. Lightly tape the tube to the patient's cheek while you're confirming placement and record the inserted distance.


- Never retract the guidewire during insertion.
- Coughing may indicate passage of the tube into the trachea. If this is suspected, remove the tube and reinsert it once the patient is comfortable.
- If resistance is encountered, remove the tube. Notify physician.
- Particular care should be taken if an endotracheal tube is in place, as it may tend to guide the feeding tube into the trachea.

Before correct placement is confirmed do not administer any fluids through the tube:


- Do not flush the tube
- Do not lubricate the internal guidewire Follow below instruction to confirm correct gastric placement

POSITION CONFIRMATION


10. Gastric aspiration

Use the syringe to obtain an aspirate of the gastric contents through the ENFit™ connector on the guidewire.

11. pH measurement


Measure the pH value of the aspirate. A $\text{pH} \leq 5.5$ indicates a correct gastric placement. Local guidelines can differ in terms of the acceptable pH range for confirming nasogastric tube placement. Always consult your local guidelines.

12. X-ray confirmation

In case no aspirate can be obtained, or the measured pH value is 5.5 or higher; Confirm gastric placement by X-ray.


△ Please follow the local radiation protection procedures.

13. Flush

Once correct positioning has been confirmed, administer at least 20ml H₂O (see 'suggested flush volume' table for volume sensitive patients) to remove gastric contents from the tube.

△ Correct placement must be confirmed before any fluid is passed through the nasogastric tube.

14. Remove the guidewire completely at this stage

△ Guidewire should never be replaced or reintroduced during or after placement or while the tube is still in the patient.

High


One o
who a

1. P


2. P

Careg
to the
enter
misce
Meas
exam

15. Close the NG-tube with the attached cap.

16. Secure the tube.

Clean the skin with alcohol to remove surface oil. The tube is gently looped and taped to patient's cheek to allow access to mouth while removing the tube from the patient's line of vision. Skin protectants could be used before taping. Attach administration set and begin feeding per physician's order and usual institution protocol.

High risk patients

One of the following procedures may be recommended for safe placement in high risk patients, e.g., those who are intubated, unconscious, have minimal or absent gag reflex, or are otherwise debilitated:

1. Placement under fluoroscopy
2. Placement under endoscopy
3. Placement under 2-stage X-ray
 - Stage 1: Confirm tube has not entered the bronchial tree.
 - Stage 2: Confirm final placement.

- Select the correct tube length. It is not allowed to alter or shorten the tube manually as the cutting edge on the tube can cause serious injuries and alter performance of the device outside of safe parameters.
- Sedation may be useful in selected patient groups; extreme care must be used as all forms of sedation have the potential to impair gag and/or swallow reflexes and increase the risk of gastric tube misplacement.

MEDICATION ADMINISTRATION AND AFTER CARE

Caregivers should apply a risk-based approach before the nasogastric route is selected to administer medication to the patient. Appropriate ancillaries including syringes should be used for medication administration through enteral feeding tubes using ENFit™ connectors of a recognized standard (ISO80369-3) in order to avoid misconnection errors.

Measures shall be taken to ensure correct drug dosing when drugs are administered through enteral tubes, for example when using (low-dose tip) ENFit™ syringes.

GENERAL FLUSHING AND SYRINGE USE POLICY

- Routinely flush the tube with water before and after the administration of feed or medication and at least every 8 hours to prevent blockage of the tube.
- As practicable, use an enteral syringe (minimum size of 20mL) for flushing the tube and always flush in a controlled manner.
- For medication administration directly into the tube, always administer using an appropriate size enteral syringe for the volume of medication but ensure delivery is in a controlled manner. In case more than one medication is administered at one time; make sure the tube is flushed thoroughly between types.
- Never use excessive force. If there is any resistance, stop the procedure and consult your clinician.

Check the position (see IFU: POSITION CONFIRMATION step 10-12) of the tube before use and each time you change the tube feed, when in doubt whether the tube is in the correct position, and at least every 8 hours.

Actual flush volumes are to be determined by local practice and locally determined clinical procedures. This table provides an indication only ensuring the tube contents (inner lumen) are completely flushed.

Tube sizes	Minimum flush volume (volume sensitive patients)
Pediatrics	0,5mL (CH4) - 2mL (CH6)
Adult sizes	20mL (CH8) - 50mL (CH14)

REMOVAL

It is recommended to (re)place the Flocare® nasogastric feeding tube every 6 weeks to minimize the risk of complications caused by material degradation or clogging within the tube. Replacement or removal of the product may differ from this recommendation, based on medical scrutiny, viability of the product and/or potential patient inconveniences'. The Flocare® nasogastric feeding tubes can be removed by gently pulling the tube out of the patient's nose. Removal should always be performed under medical supervision. Once removed from the patient the nasogastric feeding tube should be disposed in accordance with hospital, administrative and/or local government policy.

POTENTIAL COMPLICATIONS

The most commonly reported complications of nasogastric feeding are aspiration and pneumonia. Irritation around the nostrils can also occur as caused by mechanical irritation of the tube. General complications for enteral feeding are gastro-intestinal and metabolic.

SINGLE PATIENT USE

Do not reuse this product - Do not reuse, reprocess or re-sterilize this device. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

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